



How to Install the Battery

To install the battery, take these steps:

1. Hold the activity monitor as illustrated in photo #1 at left. (Make sure your index fingers are resting on the clip, not the top of the cover you're about to slide off.)
2. Apply pressure to your thumbs, pushing the cover off the NL-2000. (This is illustrated in photo #2 at left.) Set the cover aside for now.
3. Open the activity monitor.
4. Insert your right thumbnail into the groove along the right side of the device, as shown in photo #4 at left. (If you can't find the groove, see the illustration on page 17.)
5. Pull out the battery tray with your right fingers. (You may have to jiggle this a bit, or twist your fingernail slightly to get the tray open.)
6. Set the battery inside the tray, with the + sign facing upright, as illustrated in photo #6 at left.
7. Carefully slide the tray back into the NL-2000, as illustrated in photo #7 at left. Push until you hear the tray click into place. Slide the cover back on.

Note: The LCD panel will display all characters for a few seconds, and then revert to 0:00, with the first '0' flashing. The NL-2000 is now in Start-up mode. You'll need to re-enter the Start-up data (see First Things First, starting with page 19, #5).